

Secure-A-Drug Mail-Back Instructions




DO NOT THROW AWAY ANYTHING. RETAIN BOX FOR RETURN SHIPMENT.

1

Ensure fluid **solidifier** has been added to the containers and the lids are locked.

- | | |
|---------------|---------------|
| 2L container | 1 x 21g pouch |
| 4L container | 1 x 2oz pouch |
| 9L container | 3 x 2oz pouch |
| 18L container | 5 x 2oz pouch |



2

Line the box with the plastic liner. Place the Secure-A-Drug containers inside the liner and use the zip tie to seal the liner and ensure it is water tight.



3

Close the box, remove the old shipping labels and **seal the box** with three strips of shipping tape as shown below:



4

Place the adhesive return label on top of the box.

If you do not have a daily UPS pick-up service:
Schedule a UPS pick-up online at ups.com or by calling 800-742-5877



Things to Remember

Each Secure-A-Drug mailback kit includes:

- Adhesive return shipping label (if shipped via UPS)
- Plastic liner and zip tie
- 3 strips of pre-cut shipping tape
- Solidifier pouches

To re-order:

- Visit www.secureadrug.com or
- Email sales@secureadrug.com

If a spill occurs during shipping, please call Chemtrek at: 1 800 424 9300